

SUBSCRIBE:



e-newsletter



magazines

About Us | Site Map

advanced

RDH

Volume 29 Issue 6
June, 2009

- Current Issue
- Issue Archive

Home > Issue Archive > Dental Economics > Article Display



Print this article



Email this article

Text Size



BOOKMARK



Add RSS Feed

Six superior selections

March 15, 2005

Joe Blaes, DDS

FIT CHECKER™ II NDS by GC AMERICA: Finally, a new, more convenient material for checking high spots and pressure points in inlays, crowns, bridges, and dentures. FIT CHECKER II is a white-colored, vinyl polysiloxane material in an automix cartridge that will show you high spots that need to be relieved. It has plenty of working time for dentures, with a better flow under pressure. It will easily peel away from tooth and tissue surfaces. One caution - FIT CHECKER II is not compatible with the standard FIT CHECKER and should not be used on the surface of silicone-based soft relines. Order from your dealer or visit www.gcamerica.com for more information.

NuCare™ Root Conditioner by Sunstar BUTLER: Here's a new word for you - biomimetic - the study of the formation, structure, or function of biologically produced substances and materials and biological mechanisms and processes especially for the purpose of synthesizing similar products by artificial mechanisms which mimic natural ones. NuCare root conditioner uses a process that mimics the natural process to yield a natural result - biomimetic tubule occlusion to desensitize tooth root dentin. When the material comes in contact with saliva or water, it reacts and releases calcium and phosphorous ions. These ions immediately precipitate back onto the surface of the tooth to form a hydroxy-carbonate-apatite layer (HCA) which is chemically the same as the existing tooth mineral. It is the only desensitizing product cleared for subgingival use. Now you can finally do something about postoperative scaling and root-planing sensitivity. NuCare is easy to use and will make your patients happy! Order from your dealer or from the Web site at www.jbutler.com.

Quad-Tray Xtreme™ by CLINICIAN'S CHOICE Dental Products: A great product in the past has now been improved to make it even better. The distal bar that prevents the tray from flexing has been thinned by 25 percent, eliminating retro molar pad impingement and tray distortion. The lingual wall of the tray is shorter. New upper and lower retention bars lock the impression material to the tray. The mesh has been attached and extended to ensure an accurate impression of up to four teeth. You will absolutely love the new Quad-Tray Xtreme! To order, call (800) 265-3444 or visit the Web site at www.clinicianschoice.com.

Digital Dentistry by Dr. Thomas K. Hedge: This comprehensive guide to digital photography has been completely updated by one of the masters. Dr. Tom Hedge will guide you through digital photography and its applications in the dental office. Are you confused by which camera to buy? Tom will give you the information you need to make an informed choice. Also available is a terrific Digital Dentistry CD that includes PowerPoint presentations for patients as well as many forms that Tom uses in

NEW!

Ultrasonic scaling just got three times easier.

1 Dual LED Illumination 2 Ergonomic Swivel 3 Optimized Water Flow



insight
ULTRASONIC INSERTS

[Learn More](#)

Dental Economics

Looking for more news and information? Search our archives. [Click here!](#)

RECENTLY ARCHIVED ARTICLES

867-5309

Occlusion, Function, and Parafunction: Understanding the Dynamics of a Healthy Stomatognathic System
A boom for general dentists: The paradigm has shifted
A breakthrough in implant dentistry
Answering common endodontic questions: Part 2

Interested in a subscription to Dental Economics Magazine? [Click Here](#) to subscribe!

[XML](#) [RSS](#) [Click Here](#) to subscribe to the Dental Economics Feed.

RDH Magazine introduces



RDH View

a free online sampling program just for you, the Hygienist!

FREE PRODUCTS

Resource Center

- Oral Cancer Prevention Center
- Dental Continuing Education
- Dental Practice Management
- Under One Roof Blog
- From The Editors Blog
- Download Center
- Buyers Guides
- Media Kit
- White Papers
- MSDS Sheets
- Application / Technique Guides
- Classified Ads
- Free Samples

RDH Magazine introduces



RDH View

a free online sampling program just for you, the Hygienist!

FREE PRODUCTS

Did you know that acid erosion is a growing problem in children?

Introducing Sensodyne® Pronamel™ for Children
Gives kids a fighting chance against acid erosion

Online B.S. completion for licensed hygienists

- keep working
- live at home
- relocate without worry

[Click to learn more](#)



his office to be more productive. I would advise getting both! To order your copy, call (513) 777-7017 or e-mail corine@dentalhealthcenter.com.

CPR™ from IC Care:

Cosmetic Polishing Restorative™ is a unit-dose polishing agent which can be applied with a prophy brush to renew composite, porcelain, veneers, and other cosmetic restorations. We have found this product to be a great help for our hygienists, who are naturally worried about what they use to polish cosmetic restorative work. We use it as a final step in seating restorations to bring up the final luster and to minimize future staining, and the hygienist uses it to polish all cosmetic restorations. CPR is inexpensive, the only one that we found in a unit dose, and one of the few polishing agents that can be used for ceramics and composites. To order, call your dealer or visit www.iccare.net for more information.

DE

ON THE WEB

The air is salty, the islands deserted, the food delicious and - best of all - the cell phones absent. This is life on the American Eagle, an authentic Maine windjammer and National Historic Landmark. Seventy-five years ago, the ship was hauling fish off the shores of Gloucester. Today, having been lovingly and expertly refurbished by owner/captain John Foss, the ship carries passengers on four- to six-day cruises around Penobscot Bay, located in the midcoast of Maine about 90 miles north of Portland. It is, as they say in the Northeast, "wicked nice." Learn more in February on the Dental Economics Web site.

[Add RSS Feed](#)

[More Articles >](#)

[Search Archives >](#)

PennWell Dental Group Article Categories:

[Dental News](#)

[Dental Headlines](#)

[Video News](#)

[Dental Products](#)

[Industry News](#)

[Current Issue Table of Contents](#)

[Search Products Buyer's Guide >](#)

[Magazine & E-Newsletter Subscriptions >](#)



RDH eVillage®

Join the RDH community online

[PennWell Dental Group](#) | [Dental Economics](#) | [RDH](#) | [Woman Dentist Journal](#) | [Dental Office](#) | [Grand Rounds](#) | [Proofs](#) | [About Us](#) | [Contact Us](#) | [Media Kit](#)

[Dental Practice Management](#) | [Dental Continuing Education](#) | [Dental Events](#) | [PennWell](#) | [PennWell Sites](#) | [PennWell Events](#) | [Courage and Valor](#)

[Privacy Policy](#) | [Terms & Conditions](#) | [Site Map](#) | [Webmaster](#) | [RSS Feed](#) [XML](#) [RSS](#)

Copyright © 2007: PennWell Corporation, Tulsa, OK: All Rights Reserved.